

Learn more about EasyOne Air at nddmed.com 300 Brickstone Square, Suite 604 Andover, MA 01810 | (978) 470-0923 | customerservice@nddmed.com



Simple. Certain. Proven.





The original ultrasonic flow measurement is highly accurate in all flow ranges, independent of gas composition, pressure, temperature and humidity. ndd TrueFlow is a hygienic and resistance-free solution that does not require calibration during its lifetime.



EasyOne Connect software offers one database and one platform for all EasyOne point-of-care-solutions. Data management has never been easier.

Specifications

Test Types FVC, FVL, SVC, MVV and Pre/Post Testing

Languages English, French, German, Spanish

Printing Options Direct to printer or with EasyOne Connect software

EMR connectivity EasyOne Connect

Dimensions 3.4x6.1x1.4", 13 oz

Meets ATS/ERS 2005, OSHA and SSD requirements

Part Numbers

2500-2A EasyOne Air spirometry system

5050-50 EasyOne Flow Tube box of 50

5050-200 EasyOne Flow Tube box of 200